



HTA HT3000A Liquid Autosampler

Agilent 7683 Liquid Autosampler

About

The NEW HTA HT3000A Liquid Autosampler is easy to install, easy support, doesn't require additional bench space and is offered at a great price! It can be configured to serve up to two injectors and the rotating tower leaves the injector port free for manual injections or maintenance. Because of its flexibility, the HT3200A can easily be installed on all the GC and GC/MS systems available on the market!

The Agilent 7683 series automatic liquid sampler has adequate performance and reliability. This generation Agilent automatic liquid sampler is only compatible with Agilent GC Series. This is the automatic sampler to have for: accuracy and ease of use; superior reproducibility, with minimal rework; easy accommodation to changing needs; injection flexibility; and enduring reliability.

Sample Size

209 2ml Samples (2 trays)

100 Sample Tray

Options

Liquid, Headspace

Liquid Only

Injection Time

100ms

0.1 s (Standard) 4s (Cool on-column)

Requires Tray Installation

No

Yes

Benchtop Space Requirements

Does not require additional bench space

Tray requires more room on GC

Injection Ports

2

1

Condition

New

Refurbished

Multiple Injection

Yes

Yes

Chemstation Control

Yes

Yes

Warranty

2 Year Warranty

90 Day